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Silikonové endokardiální elektrody Riata a Riata ST společnosti St. Jude Medical, modely 1560, 1561, 1562, 1570, 1571, 1572, 1580, 1581, 1582, 1590, 1591, 1592, 7000, 7001, 7002, 7010, 7011, 7040, 7041, 7042

15. prosince 2010

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SEZNAM LITERATURY

- ¹ Mehta D, Nayak H, Singson M, Chao S, Jorge E, Camunas L, Gomes J. Late Complications in Patients with Pectoral Defibrillator Implants with Transvenous Defibrillator Lead Systems: High Incidence of Insulation Breakdown. *PACE*, 1998; 21:1893-1900.
- ² De Lurgio D, Sathavorn C, Mera F, Leon A, Walter P, Langberg J. Incidence and Implications of Abrasion of Implantable Cardioverter-Defibrillator Leads. *American Journal of Cardiology*, 1997, 79:1409-1411.
- ³ Pavia S, Saliba W, Wilkoff B. Lead System Dysfunction, Diagnosis, and Therapy. In: Pacifico A. ed-in-chief. *Implantable Defibrillator Therapy: A Clinical Guide*. Norwell, Ma. Kluwer Academic Publishers 2002:259-278.
- ⁴ Kleemann T, Becker T, Doenges K, Vater M, Senges J, Schneider S, Saggau W, Weisse U, Seidl K. Annual Rate of Transvenous Defibrillation Lead Defects in Implantable Cardioverter-Defibrillators over a Period > 10 years. *Circulation*. 2007;115:2474-2480.
- ⁵ Blommaert D, Louagie Y, Eucher P, Collet B, Mancini I, Govarets G, Leclercq C, De Roy L. Are Cardioverter-Defibrillator Leads Reliable Over time? *NASPE Abstracts*, 2003, 26:1123.
- ⁶ Degeratu F, Khalighi K, Peters R, Shorofsky S, Gold M. Sensing Lead Failure in Implantable Defibrillators: A Comparison of Two Commonly Used Leads. *J Cardiovasc Electrophysiol*, 2000, 11:21-24.
- ⁷ Valk S, Luijten R, Jordaens L. Insulation Damage in a Shock Wire: An Unexpected Fluoroscopic Image. *PACE*. 2010; 1-3.
- ⁸ Richards M, Warren C, Anderson M. Late failure of a single-coil Transvenous implantable cardioverter-defibrillator lead associated with conductor separation. *Europace Advanced Access* published March 27, 2010.
- ⁹ Jalal Z, Derval N, Ploux S, Bordachar P. Unusual failure of a multilumen, small-diameter implantable cardioverter-defibrillator lead. *Heart Rhythm*, 2009, Articles in Press.
- ¹⁰ Duray G, Israel C, Schmitt J, Hohnloser S. Implantable cardioverter-defibrillator lead disintegration at that level of the tricuspid valve. *Heart Rhythm Society*, 2008.
- ¹¹ HRS/EHRA Expert Consensus on Monitoring Cardiovascular Implantable Electronic Devices (CIED) dated April 2008
- ¹² Wilkoff et al., Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management, *Heart Rhythm*, Vol 6, No 7, July 2009.
- ¹³ Maisel et al., Recommendations from the Heart Rhythm Society Task Force on Lead Performance Policies and Guidelines, *Heart Rhythm*, Vol 6, No 6, June 2009.