Askina® Foam

cs	Pěnový obvaz
de	Schaumstoff-Wundauflage
en	Foam dressing
es	Vendaje de espuma
et	Ääriseva vahtkate
fi	Vaahtoside
fr	Pansement en mousse
it	Medicazione in schiuma
nl	Foamverband
no	Skumbandasje
pt	Penso de espuma
ro	Pansament din spumă
sk	Penové krytie
sv	Skumförband
th	ິ່ສົດປິດແພລໜີນິດໄວມ

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Askina® Foam

Winner Medical Co. Ltd.
Winner Industrial Park,
No 660 Bulong Road,
Longhua District,
518109 Shenzhen, China

EC REP
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80,
20537 Hamburg, Germany
B. Braun Avitum AG
Schwarzenberger Weg 73-79
34212 Melsungen, Germany

08/22

INSTRUCTIONS FOR USE

1. PRODUCT DESCRIPTION

Askina® Foam, an absorbent dressing that can reduce pain and trauma during dressing changes. The dressing absorbs wound exudate through the polyurethane foam and provides a moist wound healing environment to promote wound healing by the presence of a vapour permeable, water and bacteria resistant polyurethane film outer layer.

2. INDICATIONS

Askina® Foam dressing is designed for a wide range of moderate to highly exuding wounds such as pressure ulcers, leg and foot ulcers, and traumatic wounds e.g. first and second degree burns, and surgical wounds.

3. FUNCTION PRINCIPLES

Askina® Foam consists of polyurethane film and polyurethane foam. The dressing effectively absorbs wound exudate through the polyurethane foam and provides a moist wound healing environment to promote wound healing by the presence of a vapour permeable, water and bacteria resistant polyurethane film outer layer.

4. CONTRAINDICATIONS

- 1) Do not use as a primary treatment on 3rd degree burns.
- 2) Do not use on the wound if there are signs of infection and consult a health care professional.
- 3) Do not use on patients with a known hypersensitivity to the product or to its components.

5. INSTRUCTION FOR USE

1. Clean the wound area and instruments in accordance with normal procedures. Ensure the peri-wound skin is dry.

2. Select a size of Askina® Foam dressing enough to overlap the wound edges by at least 2 - 3 cm. if necessary, cut to shape.

3. Gently apply the foam side (white side or yellowish side) directly to the wound surface and fix it.

4. Askina® Foam may be secured by any secondary dressing (clinician preference) e.g. roll gauze, elastic or cohesive wraps, film dressing, tubular soft cloth surgical tape.

Note: The dressing change interval may be several days. Change the dressing before it is fully saturated, at signs of leakage, or as indicated by clinical practice. The same product cannot be used for more than 7 days, continuous use cannot be more than 30 days.

6. WARNINGS AND PRECAUTIONS

- 1) Do not reuse. Reuse will cause cross-contamination.
- 2) Do not use if package is damaged or open.
- 3) Do not use together with oxidizing agents such as hypochlorite solutions or hydrogen peroxide.
- 4) The product can be used by or under the supervision of health care professionals (Physicians and nurses).

7. STORAGE AND TRANSPORT CONDITIONS

1. Askina® Foam should be stored in warehouse under normal room condition, temperature and humidity of room will not be specified strictly but require cool and dry conditions. Direct sunlight shall be avoided.

2. The transportation should be keep away from direct sunlight and keep dry.

8. DISPOSAL

Disposal should be handled according to local environment procedures.

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